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unicca®: new surface implants to accelerate the osseointegration

ABSTRACT

The surface of the dental implants is a fundamental parameter to consider for improving bone regeneration based on implants. Together with the surgical techniques and materials and tools used, surface modifications of the implants have evolved over time to allow to shorten treatment times and address increasingly complex reconstructions.

The unicCa® surface is the result of the addition to the optima® surface of a layer of calcium ions. The modifications made in the development of this new surface make this advantages such as: impairment of the aging of the titanium oxide, improvement of the union bone-implant and acceleration of the phase of the osseointegration. To demonstrate the properties of this new area studies were performed in rabbit femoral condyle and in tibia of sheep. In both cases showed an increase and acceleration of the osseointegration.

The calcium, present on the the unicCa® surface ensures the cellular stimulation from the first moments after implantation until the consolidation of the tissues and the formation of the osseointegration of calcified layer which is the main constituent. This implies a regeneration peri-implant faster and of better quality.

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1. Gorbet MB, Sefton MV. Biomaterial-associated thrombosis: roles of coagulation factors, complement, platelets and leukocytes. *Biomaterials* 2004 ;25: 5681–5703.
2. Mann S. Organic matrix-mediated biomineralization. In: Mann S, editor. *Biomineralization. Principles and concepts in bioinorganic materials chemistry*. Oxford University Press, New York, 2001. 89-124.
3. Marie PJ. The calcium-sensing receptor in bone cells: a potential therapeutic target in osteoporosis. *Bone* 2010; 46: 571–576.
4. Woulfe DS. Platelet G protein-coupled receptors in hemostasis and thrombosis. *J Thromb Haemost* 2005; 3: 2193-200.
5. Gupta S, Reviakine I. Platelet activation profiles on TiO₂: effect of Ca²⁺ binding to the surface. *Biointerphases* 2012; 7: 28-40.
6. Tejero R, Rossbach P, Keller B, Anitua E, Reviakine I. Time-of-Flight Secondary Ion Mass Spectrometry with principal component analysis of titania-blood plasma interfaces. *Langmuir* 2013; 29: 902-912.
7. Anitua E, Prado R, Orive G, Tejero R. Effects of calcium-modified titanium implant surfaces on platelet activation, clot formation, and osseointegration. *J Biomed Mater Res A* 2015; 103: 969-80.
8. Sánchez-Ilárduya MB, Trouche E, Tejero R, Orive G, Reviakine I, Anitua E. Time-dependent release of growth factors from implant surfaces treated with plasma rich in growth factors. *J Biomed Mater Res A* 2013; 101: 1478-1488.
9. Gruber R, Karreth F, Kandler B, Fuerst G, Rot A, Fischer MB, Watzek G. Platelet-released supernatants increase migration and proliferation, and decrease osteogenic differentiation of bone marrow-derived mesenchymal progenitor cells under in vitro conditions. *Platelets* 2004; 15: 29-35.
10. Anitua E, Tejero R, Zalduendo MM, Orive G. Plasma Rich in Growth Factors (PRGF-Endoret) promotes bone tissue regeneration by stimulating proliferation, migration and autocrine secretion on primary human osteoblasts. *J Periodontol* 2013; 84: 1180-1190.
11. Adams CS, Mansfield K, Perlot RL, Shapiro IM. Matrix regulation of skeletal cell apoptosis. Role of calcium and phosphate ions. *J Biol Chem* 2001; 276: 20316–22.
12. Dvorak MM, Siddiqua A, Ward DT, Carter DH, Dallas SL, Nemeth EF, Riccardi D. Physiological changes in extracellular calcium concentration directly control osteoblast function in the absence of calciotropic hormones. *Proc Natl Acad Sci* 2004; 101: 5140–5.
13. Nakamura S, Matsumoto T, Sasaki J-I, Egusa H, Lee KY, Nakano T, Sohmura T, Nakahira A. Effect of calcium ion concentrations on osteogenic differentiation and hematopoietic stem cell niche-related protein expression in osteoblasts. *Tissue Eng Part A* 2010; 16: 2467–73.
14. McKee MD, Nanci A. Osteopontin at mineralized tissue interfaces in bone, teeth, and osseointegrated implants: ultrastructural distribution and implications for mineralized tissue formation, turnover, and repair. *Microsc Res Tech* 1996; 33: 141–64.
15. Davies JE, Baldan N. Scanning electron microscopy of the bone-bioactive implant interface. *J Biomed Mater Res* 1997; 36: 429–40.
16. Tengvall P. Proteins at Titanium Interfaces. In: Brunette DM, Tengvall P, Textor M, Thomsen P, editors. *Titanium in medicine*. Berlin: Springer-Verlag; 2001, 458–83.
17. Ellingsen JE. A study on the mechanism of protein adsorption to TiO₂. *Biomaterials* 1991; 12: 593–596.
18. Aita H, Hori N, Takeuchi M, Suzuki T, Yamada M, Anpo M, Ogawa T. The effect of ultraviolet functionalization of titanium on integration with bone. *Biomaterials* 2009 ;30: 1015-25.
19. Ueno T, Yamada M, Suzuki T, Minamikawa H, Sato N, Hori N, Takeuchi K, Hattori M, Ogawa T. Enhancement of bone-titanium integration profile with UV-photofunctionalized titanium in a gap healing model. *Biomaterials* 2010 ;31:1546-57.
20. Yamada Y, Yamada M, Ueda T, Sakurai K. Reduction of biofilm formation on titanium surface with ultraviolet-C pre-irradiation. *J Biomater Appl* 2013; 23: 161-171.
21. Gallardo-Moreno AM, Pacha-Olivenza MA, Fernández-Calderón MC, Pérez-Giraldo C, Bruque JM, González-Martín ML. Bactericidal behaviour of Ti6Al4V surfaces after exposure to UV-C light. *Biomaterials* 2010; 31: 5159–68.
22. Tejero R, Bierbaum S, Douglas T, Reinstorf A, Worch H, Scharnweber D. Glucuronic acid and phosphoserine act as mineralization mediators of collagen I based biomimetic substrates. *J Mater Sci Mater Med* 2010; 21: 407-18.
23. Fowler BO. Infrared studies of apatites. I. Vibrational assignments for calcium, strontium, and barium hydroxyapatites utilizing isotopic substitution. *Inorg Chem* 1974; 13: 194–207. 56.
24. Gadaleta SJ, Paschalis EP, Camacho NP, Betts F, Mendelshon R, Boskey AL. Fourier Transform infrared spectroscopy of synthetic and biological apatites. In: Amjad Z, editor. *Mineral scale formation and inhibition*. New York: Plenum Press; 1995; 57:283–294.
25. Klee WE, Engel G. Infrared spectra of the phosphate ions in various apatites. *J Inorg Nucl Chem* 1970; 32: 1837–1843.
26. Anitua E. Perimplantitis: Un nuevo enfoque a la prevención y su tratamiento. Ed. Eduardo Anitua. Teamwork Media, Vitoria, España, 2014.
27. WP Superficie optima®: Balance entre osteointegración y riesgo de peri-implantitis. Estudio celular. Teamwork Media, Vitoria, España, 2013.
28. Anitua E, Piñas L, Murias A, Prado R, Tejero R. Effects of calcium ions on titanium surfaces for bone regeneration. *Colloids Surf B Biointerphases* 2015 ;130: 173-81.
29. Anitua E, Prado R, Orive G, Tejero R. Effects of calcium-modified titanium implant surfaces on platelet activation, clot formation, and osseointegration. *J Biomed Mater Res A* 2015; 103: 969-80.
30. Anitua E, Tejero R, Pacha-Olivenza MÁ, Fernández-Calderón MC, Delgado-Rastrolo M, Zalduendo MM, Troya M, Pérez-Giraldo C, González-Martín ML. Balancing microbial and mammalian cell functions on calcium ion-modified implant surfaces. *J Biomed Res B Appl Biomater* 2018; 106 (1): 421-432.
31. Anitua E, Murias-Freijo A, Piñas L, Tejero R, Prado R, Orive G. Nontraumatic Implant Explantation: A Biomechanical and Biological Analysis in Sheep Tibia. *J Oral Implantol* 2016; 42: 3-11.